



defrost tube racks in less than 20 minutes

The choice for academic, biotech and research labs to reduce time intensive thawing processes.

fast performance

- Defrosts 96 tubes in under 20 minutes
- Standalone system
- No warm up time


small footprint

- Saves benchtop space

best price

- under \$650

Thawing Station



specifications	electrical	dimensions	weight
	110-220 VAC 50/60HZ	6.25"w x 4.25"d x 11"h 16 cm x 11 cm x 28 cm	5 lbs. (2.3 kg)
Compatible Labware: ANSI/SLAS tube racks from major manufacturers including Biosigma, FluidX, LVL, Micronic, Thermo Scientific (ABgene, Matrix, Nunc) and others			

BioMicroLab

sptlabtech.com

sptlabtech